

HDBS0637

Product Information

Product SKU:	HDBS0637	Host Species:	Whole mAb	Size:	100 µg
Target:	PD-1	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	Q15116
Isotype:	IgG4		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1

Common Name: BMS-936558;BMS936558;BMS 936558;MDX-1106;MDX1106;MDX 1106;ONO-4538;ONO4538;ONO 4538;OPDIVO;OPDUALAG (combination of Nivolumab and Relatlimab)

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.