

**HDBS0646**

---

## Product Information

<b>Product SKU:</b> HDBS0646	<b>Host Species:</b> Whole mAb	<b>Size:</b> 100 µg
<b>Target:</b> VEGFR2	<b>Reactivity:</b> Human	

---

## Additional Information

<b>Conjugate:</b> Unconjugated	<b>Uniprot:</b> P35968
<b>Isotype:</b> IgG1	

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** FLK1;CD309;VEGFR;VEGFR2

**Common Name:** SSS-23;SSS23;SSS 23;TTAC-0001;TTAC0001;TTAC 0001;TTAC-1;TTAC1;TTAC 1

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.