

HDBS0655

Product Information

Product SKU: HDBS0655	Host Species: Whole mAb	Size: 100 µg
Target: IL6	Reactivity: Human	

Additional Information

Conjugate: Unconjugated	Uniprot: P05231
Isotype: IgG2	

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CDF;HGF;HSF;BSF2;IL-6;BSF-2;IFNB2;IFN-beta-2

Common Name: PF-04236921;PF-4236921;PF 04236921;PF 4236921;PF04236921;PF4236921;TOUR-006;TOUR006;TOUR 006

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.