

HDBS0658

Product Information

Product SKU: HDBS0658 **Host Species:** Whole mAb **Size:** 100 µg
Target: TLR4 **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** O00206
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: TOLL;CD284;TLR-4;ARMD10

Common Name: EB 05;EB05;EB-05;EB 5;EB-5;EB5;EB-07;EB07;EB-07;EB-7;EB7;EB 7;NI-0101;NI 0101;NI0101;NI-101;NI101;NI 101

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.