

HDBS0662

Product Information

| | | | | | |
|---------------------|----------|----------------------|-----------|--------------|--------|
| Product SKU: | HDBS0662 | Host Species: | Whole mAb | Size: | 100 µg |
| Target: | PD-1 | Reactivity: | Human | | |

Additional Information

| | | | |
|-------------------|--------------|-----------------|--------|
| Conjugate: | Unconjugated | Uniprot: | Q15116 |
| Isotype: | IgG1 | | |

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1

Common Name: LY-3462817;LY3462817;LY 3462817

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.