

HDBS0663

Product Information

Product SKU: HDBS0663 **Host Species:** Whole mAb **Size:** 100 µg
Target: Her2 **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** P04626
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: NEU;NGL;HER2;TKR1;CD340;HER-2;VSCN2;MLN 19;MLN-19;c-ERB2;c-ERB-2;HER-2/neu;p185(erbB2)

Common Name: R-1273;R1273;R 1273;RG-1273;RG 1273;RG1273;RO-4368451;RO4368451;RO 4368451

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.