

HDBS0668

Product Information

Product SKU: HDBS0668	Host Species: Whole mAb	Size: 100 µg
Target: CD123	Reactivity: Human	

Additional Information

Conjugate: Unconjugated	Uniprot: P26951
Isotype: IgG1	

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: IL3R;CD123;IL3RX;IL3RY;IL3RAY;hIL-3Ra;IL-3R-alpha

Common Name: IMGN-632;IMGN632;IMGN 632

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.