

HDBS0673

Product Information

Product SKU: HDBS0673	Host Species: Whole mAb	Size: 100 µg
Target: SNCA	Reactivity: Human	

Additional Information

Conjugate: Unconjugated	Uniprot: P37840
Isotype: IgG1	

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: PD1;NACP;PARK1;PARK4

Common Name: NEOD-002;NEOD002;NEOD 002;NEOD-2;NEOD2;NEOD 2;PRX-002;PRX002;PRX 002;PRX-2;PRX2;PRX 2;RG 7935;RG7935;RG-7935;RO 7046015;RO-7046015;RO7046015

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.