

HDBS0693

## Product Information

**Product SKU:** HDBS0693      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** PD-1      **Reactivity:** Human

## Additional Information

**Conjugate:** Unconjugated      **Uniprot:** Q15116  
**Isotype:** IgG4

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1

**Common Name:** AEX-1188;AEX1188;AEX 1188;INCMGA-0012;INCMGA 0012;INCMGA0012;INCMGA12;INCMGA-12;INCMGA 12;MGA-012;MGA012;MGA 012;MGA-12;MGA12;MGA 12;Zynyz

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.