

**HDBS0713**

---

## Product Information

**Product SKU:** HDBS0713      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** SOST      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Unconjugated      **Uniprot:** Q9BQB4  
**Isotype:** IgG2

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CDD;VBCH;DAND6;SOST1

**Common Name:** BPS-804;BPS804;BPS 804;MOR-05813;MOR05813;MOR 05813;MOR-5813;MOR5813;MOR 5813

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.