

HDBS0717

Product Information

Product SKU: HDBS0717 **Host Species:** Whole mAb **Size:** 100 µg
Target: PD-L1 **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** Q9NZQ7
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1

Common Name: STI-1014;STI1014;STI 1014;STI-A1014;STI A1014;ZKAB-001;ZKAB001;ZKAB 001;ZKAB-1;ZKAB1;ZKAB 1

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.