

**HDBS0721**

---

## Product Information

<b>Product SKU:</b>	HDBS0721	<b>Host Species:</b>	Whole mAb	<b>Size:</b>	100 µg
<b>Target:</b>	CD40	<b>Reactivity:</b>	Human		

---

## Additional Information

<b>Conjugate:</b>	Unconjugated	<b>Uniprot:</b>	P25942
<b>Isotype:</b>	IgG1		

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CD40;p50;Bp50;CDW40;TNFRSF5

**Common Name:** APX 005;APX005;APX 005;APX5;APX-5;APX 5;EPI-0050;EPI0050;EPI 0050;EPI 50;EPI-50;EPI50

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.