

HDBS0734

Product Information

| | | | | | |
|---------------------|----------|----------------------|-----------|--------------|--------|
| Product SKU: | HDBS0734 | Host Species: | Whole mAb | Size: | 100 µg |
| Target: | NGF | Reactivity: | Human | | |

Additional Information

| | | | |
|-------------------|--------------|-----------------|--------|
| Conjugate: | Unconjugated | Uniprot: | P01138 |
| Isotype: | IgG2 | | |

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: NGFB;HSAN5;Beta-NGF;NGF;NGFB;-NGF

Common Name: PF-04383119;PF04383119;PF-4383119;PF4383119;RI-624;RI624;RI 624;RN-624;RN624;RN 624

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.