

HDBS0750

Product Information

Product SKU: HDBS0750 **Host Species:** Whole mAb **Size:** 100 µg
Target: PD-1 **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** Q15116
Isotype: IgG4

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1

Common Name: CHS-007;CHS007;CHS 007;CHS-7;CHS7;CHS 7;JS 001;JS001;JS-001;JS1;JS-1;JS 1;SO-001;SO 001;SO001;SO1;SO-1;SO 1;TAB-001;TAB001;TAB 001;TAB1;TAB-1;TAB 1;Tripleitriumab;Triprizumab;Tuoyi

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.