

HDBS0751

Product Information

Product SKU: HDBS0751	Host Species: Whole mAb	Size: 100 µg
Target: IL33	Reactivity: Human	

Additional Information

Conjugate: Unconjugated	Uniprot: O95760
Isotype: IgG4	

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: DVS27;IL1F11;NF-HEV;NFEHEV;C9orf26

Common Name: LY-3375880;LY3375880;LY 3375880;ZB-880;ZB880;ZB 880

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.