

HDBS0754

---

## Product Information

**Product SKU:** HDBS0754      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** CTLA-4      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Unconjugated      **Uniprot:** P16410  
**Isotype:** IgG2

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CD;GSE;GRD4;ALPS5;CD152;CTLA-4;IDDM12;CELIAC3

**Common Name:** CP-675;CP675;CP 675;CP-675206;CP 675206;CP675206;IMJUDO;Ticilimumab

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.