

HDBS0762

Product Information

Product SKU: HDBS0762 **Host Species:** Whole mAb **Size:** 100 µg
Target: HGF **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** P14210
Isotype: IgG4

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: SF;HGFB;HPTA;F-TCF;DFNB39

Common Name: YYB101;YYB-101;YYB 101;CLM-101;CLM 101;CLM101;NOV 1105;NOV-1105;NOV1105

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.