

**HDBS0786**

---

## Product Information

|                     |          |                      |           |              |        |
|---------------------|----------|----------------------|-----------|--------------|--------|
| <b>Product SKU:</b> | HDBS0786 | <b>Host Species:</b> | Whole mAb | <b>Size:</b> | 100 µg |
| <b>Target:</b>      | APRIL    | <b>Reactivity:</b>   | Human     |              |        |

---

## Additional Information

|                   |              |                 |        |
|-------------------|--------------|-----------------|--------|
| <b>Conjugate:</b> | Unconjugated | <b>Uniprot:</b> | O75888 |
| <b>Isotype:</b>   | IgG4         |                 |        |

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** APRIL;CD256;TALL2;ZTNF2;TALL-2;TNLG7B;TRDL-1;UNQ383/PRO715

**Common Name:** BION-1301;BION1301;BION 1301

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.