

### Product Information

<b>Product SKU:</b>	HDBS0788	<b>Host Species:</b>	Whole mAb	<b>Size:</b>	100 µg
<b>Target:</b>	PD-1	<b>Reactivity:</b>	Human		

### Additional Information

<b>Conjugate:</b>	Unconjugated	<b>Uniprot:</b>	Q15116
<b>Isotype:</b>	IgG4		

### Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1

**Common Name:** AB 122;AB-122;AB122;GLS 010;GLS010;GLS010;GLS10;GLS-10;GSL 10;GS-0122;GS0122;GS 0122;GS-122;GS122;GS 122;WBP 3055;WBP3055;WBP-3055;YUTUO

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.