

HDBS0790

Product Information

Product SKU: HDBS0790 **Host Species:** Whole mAb (Canine) **Size:** 100 µg
Target: Canine IL31RA **Reactivity:** Canine

Additional Information

Conjugate: Unconjugated **Uniprot:** XP_013963900.1
Isotype: IgG2

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384; zcytoR17

Common Name: VTQ-2101;VTQ2101

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.