

HDBS0801

Product Information

Product SKU: HDBS0801 **Host Species:** Whole mAb **Size:** 100 µg
Target: INSR **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** P06213
Isotype: IgG2

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: HHF5;CD220

Common Name: RZ-358;X 358;XOMA 247;XOMA-247;XOMA247;XOMA 358;XOMA-358;XOMA358;XPA.15.247

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.