

HDBS0803

Product Information

Product SKU:	HDBS0803	Host Species:	Whole mAb	Size:	100 µg
Target:	KIR3DL2	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	P43630
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: 3DL2;p140;NKAT4;CD158K;NKAT-4;NKAT4B;KIR-3DL2

Common Name: IPH 41;IPH-41;IPH41;IPH4102;IPH-4102;IPH 4102

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.