

HDBS0810

Product Information

Product SKU: HDBS0810 **Host Species:** Whole mAb **Size:** 100 µg
Target: IL36RN **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** Q9UBH0
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: FIL1;FIL1D;IL1F5;IL1L1;PSORP;IL1HY1;IL1RP3;IL36RA;IL-36Ra;PSORS14;FIL1(Delta)

Common Name: BI 655130;BI655130;BI-655130;SPEVIGO

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.