



HDBS0819

Product Information

Product SKU: HDBS0819 **Host Species:** Whole mAb **Size:** 100 µg
Target: CXCL8 **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P10145
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: IL8;NAF;GCP1;LECT;LUCT;NAP1;GCP-1;LYNAP;MDNCF;MONAP;NAP-1;SCYB8

Common Name: BMS-986253;BMS986253;HuMax-IL8;HuMax-Inflam;MDX018;MDX-018

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.