

# Biotinylated Anti-OX40 Ligand (Amlitelimab biosimilar) mAb

HDBS0828



## Product Information

**Product SKU:** HDBS0828      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** OX40 Ligand      **Reactivity:** Human

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P23510  
**Isotype:** IgG4

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** GP34;CD252;OX40L;TXGP1;CD134L;OX-40L;TNLG2B;TNFSF4;Glycoprotein  
Gp34;CD134 ligand

**Common Name:** KY-1005;KY1005;SAR-445229

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.