

Biotinylated Anti-CSF1R (Axatilimab biosimilar) mAb

HDBS0836



Product Information

Product SKU: HDBS0836 **Host Species:** Whole mAb **Size:** 100 µg
Target: CSF1R **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P07333
Isotype: IgG4

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: FMS;CSFR;FIM2;GPSC;HDLS;C-FMS;CD115;HDLS1;CSF-1R;BANDDOS;M-CSF-R

Common Name: INCA-34176;INCA34176;SNDX-6352;SNDX6352;UCB 6352;UCB6352;UCB-6352

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.