Biotinylat mAb HDBS0873	ed Anti-l	L4RA (Comek	ibart biosi	milar)	AssayGenie	
Product Inform	ation					
Product SKU: Target:	HDBS0873 IL4RA	Host Species Reactivity:	<b>s</b> : Whole mAb Human	Size	100 µg	
Additional Infor	mation					
<b>Conjugate</b> :	Biotinylated Uniprot: P24394					
lsotype:	lgG4					
Protein Informa	tion					
Background:		Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.				
Synonyms:	nonyms: CD124;IL4RA;IL-4RA					
Common Name:		BC005;BC-005;BC 005;MG-010;MG010;MG 010;MG-K10;MG K10				
Recommended Dilutions:		-				
<b>Applications</b> :		-				
Usage: Research use only						
<b>Formulation</b> :		Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as				
protectants before lyophilization. Pleas instructions of reconstitution. <b>Storage</b> : Store at -20°C to -80°C for 12 months in not intended for use within a month, al freezing and thawing). Lyophilized prote					orm. After reconstitution, if	