

HDBS0879

Product Information

Product SKU: HDBS0879 **Host Species:** Whole mAb **Size:** 100 µg
Target: C5 **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P01031
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: C5D;C5a;C5b;ECLZB;CPAMD4

Common Name: CH 7092230;CH7092230;CH-7092230;RG 6107;RG-6107;RG6107;RO 7092230;RO-7112689;RO7112689;SKY 59;SKY-59;SKY59

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.