

# Biotinylated Anti-IL2RA (Daclizumab biosimilar) mAb

HDBS0881



## Product Information

**Product SKU:** HDBS0881      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** IL2RA      **Reactivity:** Human

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P01589  
**Isotype:** IgG1

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** p55;CD25;IL2R;IMD41;TCGFR;IDDM10;IL2-RA;IL-2-RA

**Common Name:** BIIB-019;BIIB019;Dacliximab;RO 247375;RO247375;RO-247375;Zenapax;Zinbryta

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.