

Biotinylated Anti-CD40 Ligand (Frexalimab biosimilar) mAb

HDBS0932



Product Information

Product SKU: HDBS0932 **Host Species:** Whole mAb **Size:** 100 µg
Target: CD40 Ligand **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P29965
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: IGM;IMD3;TRAP;gp39;CD154;CD40L;HIGM1;T-BAM;TNFSF5;hCD40L

Common Name: INX 021;INX021;INX-021;SAR-441344;SAR441344;SAR 441344

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.