

Biotinylated Anti-IL1B (Gevokizumab biosimilar) mAb

HDBS0938



Product Information

Product SKU: HDBS0938 **Host Species:** Whole mAb **Size:** 100 µg
Target: IL1B **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P01584
Isotype: IgG2

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: IL-1;IL1F2;IL1beta;IL1-BETA

Common Name: S-78989;S78989;S 78989;VPM 087;VPM087;VPM-087;VPM87;VPM-87;XMA-005.2;XMA-0052;XMA52;XOMA-052;XOMA052;XOMA52;XOMA-52

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.