



## HDBS0952

### Product Information

**Product SKU:** HDBS0952      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** CCR8      **Reactivity:** Human

### Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P51685  
**Isotype:** IgG1

### Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CY6;TER1;CCR-8;CKRL1;CDw198;CMKBR8;GPRCY6;CMKBRL2;CC-CKR-8

**Common Name:** BMS 986340;BMS-986340;BMS986340

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.