

HDBS0992

Product Information

Product SKU: HDBS0992 **Host Species:** Whole mAb **Size:** 100 µg
Target: CD73 **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P21589
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: NT;eN;NT5;NTE;eNT;CD73;E5NT;CALJA

Common Name: COR-004;COR004;COR 004;COR-4;COR4;COR 4;CPI-006;CPI006;CPI 006;CPI-6;CPI6;CPI 6;CPX-006;CPX006;CPX 006;CPX-6;CPX6;CPX 6

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.