

HDBS1003

Product Information

Product SKU:	HDBS1003	Host Species:	Whole mAb	Size:	100 µg
Target:	EGFR	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P00533
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: ERBB;ERRP;HER1;mENA;ERBB1;NNCIS;PIG61;NISBD2

Common Name: DE-766;DE766;DE 766;KI-0501;KI0501;KI 0501;KI-501;KI 501;KI501;KI-0502;KI0502;KI 0502;KI-502;KI502;KI 502;OSAG101;OSAG-101;OSAG 101;Taixinsheng;TheraCIM;Theraloc;VECTHIX

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.