

HDBS1030

Product Information

Product SKU: HDBS1030	Host Species: Whole mAb	Size: 100 µg
Target: SEMA4D	Reactivity: Human	

Additional Information

Conjugate: Biotinylated	Uniprot: Q92854
Isotype: IgG4	

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: A8;GR3;BB18;CD100;COLL4;SEMAJ;coll-4;C9orf164;M-sema-G

Common Name: VX15/2503;VX 15;VX-15;VX 15;VX 15/2503;VX-15/2503

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.