

HDBS1038

Product Information

Product SKU:	HDBS1038	Host Species:	Whole mAb	Size:	100 µg
Target:	LILRB2	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	Q8N423
Isotype:	IgG4		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: ILT4;LIR2;CD85D;ILT-4;LIR-2;MIR10;MIR-10

Common Name: JTX-8064;JTX8064;JTX 8064

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.