

HDBS1049

Product Information

Product SKU:	HDBS1049	Host Species:	Whole mAb	Size:	100 µg
Target:	LIGHT	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	O43557
Isotype:	IgG4		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: LTg;CD258;HVEM;LIGHT;TNFSF14

Common Name: AEVI-002;AEVI002;AEVI 002;AEVI-2;AEVI2;AEVI 2;AVTX-002;AVTX002;AVTX 002;AVTX-2;AVTX2;AVTX 2;CERC 002;CERC-002;CERC002;CERC 2;CERC-2;CERC2;MDGN-002;MDGN002;MDGN 002;MDGN 2;MDGN-2;MDGN2;SAR 252067;SAR252067;SAR-252067

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.