

**HDBS1060**

---

## Product Information

**Product SKU:** HDBS1060      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** IL5      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P05113  
**Isotype:** IgG4

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** EDF;TRF;IL-5;Interleukin-5

**Common Name:** CDP-835;CDP835;CDP 835;CEP-38072;CEP38072;CEP-38072;CINQAERO;CINQAIR;Cinquil;CTx-55700;CTx55700;Ctx 55700;DCP-835;DCP835;DCP 835;SCH 55700;SCH55600;SCH-55700;TRFK-5;TRFK 5;TRFK5

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.