

Product Information

Product SKU:	HDBS1063	Host Species:	Whole mAb	Size:	100 µg
Target:	OX40	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P43489
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: OX40;ACT35;CD134;IMD16;TXGP1L;TNFRSF4

Common Name: ABBV-368;ABBV368;ABBV 368;PR-1628103;PR1628103;PR 1628103

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.