

**HDBS1065**

---

## Product Information

**Product SKU:** HDBS1065      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** C1S      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P09871  
**Isotype:** IgG4

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** EDSPD2

**Common Name:** BIVV-020;BIVV020;BIVV 020;SAR-445088;SAR445088;SAR 445088;TNT-020;TNT020;TNT 020

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.