

HDBS1069

Product Information

Product SKU: HDBS1069 **Host Species:** Whole mAb **Size:** 100 µg
Target: IL23A **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** Q9NPF7
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: P19;SGRF;IL-23;IL-23A;IL23P19

Common Name: ABBV 066;ABBV066;ABBV 066;ABBV 66;ABBV-66;ABBV66;BI-655066;BI655066;BI655066;SKYRIZ

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.