

## HDBS1083

---

### Product Information

**Product SKU:** HDBS1083      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** APRIL      **Reactivity:** Human

---

### Additional Information

**Conjugate:** Biotinylated      **Uniprot:** O75888  
**Isotype:** IgG2

---

### Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** APRIL;CD256;TALL2;ZTNF2;TALL-2;TNLG7B;TRDL-1;UNQ383/PRO715

**Common Name:** VIS-649;VIS649;VIS 649

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.