

Product Information

Product SKU:	HDBS1099	Host Species:	Whole mAb	Size:	100 µg
Target:	CD19	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P15391
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CD19;B4;CVID3;MGC12802

Common Name: ICP-B04;ICP B04;MOR 208;MOR208;MOR-208;XENP-5574;XENP5574;XENP 5574;XmAb-5574;XmAb5574;XmAb 5574

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.