

HDBS1103

Product Information

Product SKU: HDBS1103 **Host Species:** Whole mAb **Size:** 100 µg
Target: NGF **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P01138
Isotype: IgG2

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: NGFB;HSAN5;Beta-NGF;NGF;NGFB;-NGF

Common Name: PF-04383119;PF04383119;PF-4383119;PF4383119;RI-624;RI624;RI 624;RN-624;RN624;RN 624

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.