

HDBS1112

Product Information

Product SKU: HDBS1112 **Host Species:** Whole mAb **Size:** 100 µg
Target: MAPT **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P10636
Isotype: IgG4

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: TAU;FTD1;MSTD;PPND;DDPAC;MAPTL;MTBT1;MTBT2;tau-40;FTDP-17;PPP1R103;Tau-PHF6

Common Name: ABBV-8E12;ABBV8E12;ABBV 8E12

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.