

HDBS1113

Product Information

Product SKU:	HDBS1113	Host Species:	Whole mAb	Size:	100 µg
Target:	IL23A	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	Q9NPF7
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: P19;SGRF;IL-23;IL-23A;IL23P19

Common Name: Ilumetri;ILUMYA;MK-3222;MK3222;MK 3222;SCH-900222;SCH 900222;SCH900222;SUNPG 1622;SUNPG1622;SUNPG-1622;SUNPG 1623;SUNPG-1623;SUNPG1623

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.