

# Biotinylated Anti-CD20 (Ublituximab biosimilar) mAb

HDBS1129



## Product Information

**Product SKU:** HDBS1129      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** CD20      **Reactivity:** Human

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P11836  
**Isotype:** IgG1

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** B1;S7;Bp35;CD20;FMC7;CVID5;LEU-16;MS4A1

**Common Name:** BRIUMV;EMAB-6;EMAB6;EMAB 6;EMAB-603;EMAB603;EMAB 603;LFB-R603;LFB R603;TG-1101;TG1101;RG 1101;TG-20;TG20;TG 20;TGTX-1101;TGTX1101;TGTX 1101

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.