

HDBS1133

Product Information

Product SKU: HDBS1133 **Host Species:** Whole mAb **Size:** 100 µg
Target: 4-1BB **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** Q07011
Isotype: IgG4

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: ILA;4-1BB;CD137;CDw137;IMD109 ; TNFRSF9

Common Name: BMS-663513;BMS663513;BMS 663513;BMS-66513;BMS66513;BMS 66513;ONO-4481;ONO4481;ONO 4481

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.