

**HDBS1134**

---

## Product Information

**Product SKU:** HDBS1134      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** APP      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P05067  
**Isotype:** IgG1

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** AAA;AD1;PN2;ABPP;APPI;CVAP;ABETA;PN-II;preA4;CTFgamma;alpha-sAPP

**Common Name:** PR-1905069;PR 1905069;PR1905069

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.