



mAb

HDBS1169

Product Information

Product SKU: HDBS1169 **Host Species:** Whole mAb **Size:** 100 µg
Target: MUC5AC **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P98088
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: TBM;leB;MUC5;mucin

Common Name: NEO-101;NEO101;NEO-102;NEO102;NPC-1C;NPC1C

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.